

Publications Using *Gaussia* Luciferase:

- 1: Chopra A. *Gaussia princeps* luciferase . 2008 Jan 31 [updated 2008 Feb 28]. Molecular Imaging and Contrast Agent Database (MICAD) [Internet]. Bethesda (MD): National Center for Biotechnology Information (US); 2004-2010. Available from <http://www.ncbi.nlm.nih.gov/bookshelf/br.fcgi?book=micad&part=GLuc>
PubMed PMID: 20641352.
- 2: Chen F, Xu Z, Lu J, Lü X, Mu WL, Wang YJ, Liu DP, Liang CC. *Gaussia* luciferase reporter assay for assessment of gene delivery systems in vivo. *Chin Med Sci J*. 2010 Jun;25(2):95-9. PubMed PMID: 20598231.
- 3: Sharma R, Tsuchiya M, Skobe Z, Tannous BA, Bartlett JD. The acid test of fluoride: how pH modulates toxicity. *PLoS One*. 2010 May 28;5(5):e10895. PubMed PMID: 20531944; PubMed Central PMCID: PMC2878349.
- 4: Kimura T, Hiraoka K, Kasahara N, Logg CR. Optimization of enzyme-substrate pairing for bioluminescence imaging of gene transfer using *Renilla* and *Gaussia* luciferases. *J Gene Med*. 2010 Jun;12(6):528-37. PubMed PMID: 20527045.
- 5: Andreu N, Zelmer A, Fletcher T, Elkington PT, Ward TH, Ripoll J, Parish T, Bancroft GJ, Schaible U, Robertson BD, Wiles S. Optimisation of bioluminescent reporters for use with mycobacteria. *PLoS One*. 2010 May 24;5(5):e10777. PubMed PMID: 20520722; PubMed Central PMCID: PMC2875389.
- 6: Mehraein-Ghomi F, Basu HS, Church DR, Hoffmann FM, Wilding G. Androgen receptor requires JunD as a coactivator to switch on an oxidative stress generation pathway in prostate cancer cells. *Cancer Res*. 2010 Jun 1;70(11):4560-8. Epub 2010 May 11. PubMed PMID: 20460526; PubMed Central PMCID:PMC2909607.
- 7: Rathnayaka T, Tawa M, Sohya S, Yohda M, Kuroda Y. Biophysical characterization of highly active recombinant *Gaussia* luciferase expressed in *Escherichia coli*. *Biochim Biophys Acta*. 2010 Sep;1804(9):1902-7. Epub 2010 May 7. PubMed PMID: 20452471.
- 8: Ambolet-Camoit A, Bui LC, Pierre S, Chevallier A, Marchand A, Coumoul X, Garlatti M, Andreau K, Barouki R, Aggerbeck M. 2,3,7,8-tetrachlorodibenzo-p-dioxin counteracts the p53 response to a genotoxicant by upregulating expression of the metastasis marker *agr2* in the hepatocarcinoma cell line HepG2. *Toxicol Sci*. 2010 Jun;115(2):501-12. Epub 2010 Mar 18. PubMed PMID: 20299546.

- 9: Lwa TR, Tan CH, Lew QJ, Chu KL, Tan J, Lee YY, Chao SH. Identification of cellular genes critical to recombinant protein production using a Gaussia luciferase-based siRNA screening system. *J Biotechnol.* 2010 Apr 15;146(4):160-8. Epub 2010 Feb 25. PubMed PMID: 20188772.
- 10: Feng B, Tomizawa K, Michiue H, Han XJ, Miyatake S, Matsui H. Development of a bifunctional immunoliposome system for combined drug delivery and imaging in vivo. *Biomaterials.* 2010 May;31(14):4139-45. Epub 2010 Feb 10. PubMed PMID: 20149431.
- 11: Tian W, Wang G, Luo S, Dong X, Fu X, Tan W, Wu X. [Real time ex vivo detection and dynamic monitoring of in vivo expression of secreted luciferase gene injected by hydrodynamic method]. *Sheng Wu Gong Cheng Xue Bao.* 2009 Oct;25(10):1552-7. Chinese. PubMed PMID: 20112701.
- 12: Nakanishi H, Higuchi Y, Kawakami S, Yamashita F, Hashida M. piggyBac transposon-mediated long-term gene expression in mice. *Mol Ther.* 2010 Apr;18(4):707-14. Epub 2010 Jan 26. PubMed PMID: 20104210; PubMed Central PMCID: PMC2862518.
- 13: Luo ST, Tian WH, Wang G, Dong XY, Yang L, Wu XB. [Comparison of TTR and CMV promoters in vivo and in vitro via a secreted luciferase reporter system]. *Bing Du Xue Bao.* 2009 Nov;25(6):424-9. Chinese. PubMed PMID: 20077932.
- 14: Chung E, Yamashita H, Au P, Tannous BA, Fukumura D, Jain RK. Secreted Gaussia luciferase as a biomarker for monitoring tumor progression and treatment response of systemic metastases. *PLoS One.* 2009 Dec 15;4(12):e8316. PubMed PMID: 20016816; PubMed Central PMCID: PMC2789383.
- 15: Patel KG, Ng PP, Kuo CC, Levy S, Levy R, Swartz JR. Cell-free production of Gaussia princeps luciferase--antibody fragment bioconjugates for ex vivo detection of tumor cells. *Biochem Biophys Res Commun.* 2009 Dec 18;390(3):971-6. Epub 2009 Oct 21. PubMed PMID: 19852937.
- 16: Welsh JP, Patel KG, Manthiram K, Swartz JR. Multiply mutated Gaussia luciferases provide prolonged and intense bioluminescence. *Biochem Biophys Res Commun.* 2009 Nov 27;389(4):563-8. PubMed PMID: 19825431.
- 17: Badr CE, Niers JM, Tjon-Kon-Fat LA, Noske DP, Wurdinger T, Tannous BA. Real-time monitoring of nuclear factor kappaB activity in cultured cells and in animal models. *Mol Imaging.* 2009 Sep-Oct;8(5):278-90. PubMed PMID: 19796605; PubMed Central PMCID: PMC2856067.

- 18: Enjalbert B, Rachini A, Vedyappan G, Pietrella D, Spaccapelo R, Vecchiarelli A, Brown AJ, d'Enfert C. A multifunctional, synthetic *Gaussia princeps* luciferase reporter for live imaging of *Candida albicans* infections. *Infect Immun*. 2009 Nov;77(11):4847-58. Epub 2009 Aug 17. PubMed PMID: 19687206; PubMed Central PMCID: PMC2772526.
- 19: Michelini E, Cevenini L, Mezzanotte L, Roda A. Luminescent probes and visualization of bioluminescence. *Methods Mol Biol*. 2009;574:1-13. PubMed PMID: 19685295.
- 20: Alexeyev MF, Fayzulin R, Shokolenko IN, Pastukh V. A retro-lentiviral system for doxycycline-inducible gene expression and gene knockdown in cells with limited proliferative capacity. *Mol Biol Rep*. 2010 Apr;37(4):1987-91. Epub 2009 Aug 5. PubMed PMID: 19655272; PubMed Central PMCID: PMC2832089.
- 21: Maguire CA, Deliolanis NC, Pike L, Niers JM, Tjon-Kon-Fat LA, Sena-Esteves M, Tannous BA. *Gaussia* luciferase variant for high-throughput functional screening applications. *Anal Chem*. 2009 Aug 15;81(16):7102-6. PubMed PMID: 19601604; PubMed Central PMCID: PMC2846205.
- 22: Luker K, Gupta M, Luker G. Bioluminescent CXCL12 fusion protein for cellular studies of CXCR4 and CXCR7. *Biotechniques*. 2009 Jul;47(1):625-32. PubMed PMID: 19594447.
- 23: Lerche K, Hallmann A. Stable nuclear transformation of *Gonium pectorale*. *BMC Biotechnol*. 2009 Jul 10;9:64. PubMed PMID: 19591675; PubMed Central PMCID: PMC2720962.
- 24: Yang J, Richmond AJ. Monitoring NF-kappaB mediated chemokine transcription in tumorigenesis. *Methods Enzymol*. 2009;460:347-55. PubMed PMID: 19446734.
- 25: Kim SB, Sato M, Tao H. Genetically encoded bioluminescent indicators for stress hormones. *Anal Chem*. 2009 May 15;81(10):3760-8. PubMed PMID: 19388668.
- 26: Tannous BA. *Gaussia* luciferase reporter assay for monitoring biological processes in culture and in vivo. *Nat Protoc*. 2009;4(4):582-91. PubMed PMID: 19373229; PubMed Central PMCID: PMC2692611.
- 27: Serganova I, Moroz E, Vider J, Gogiberidze G, Moroz M, Pillarsetty N, Doubrovin M, Minn A, Thaler HT, Massague J, Gelovani J, Blasberg R. Multimodality imaging of TGFbeta signaling in breast cancer metastases. *FASEB J*. 2009 Aug;23(8):2662-72. Epub 2009 Mar 26. PubMed PMID: 19325038; PubMed Central PMCID:PMC2717767.

- 28: Kaihara A, Sunami A, Kurokawa J, Furukawa T. A genetically encoded bioluminescent indicator for the sodium channel activity in living cells. *J Am Chem Soc.* 2009 Apr 1;131(12):4188-9. PubMed PMID: 19271706.
- 29: Santos EB, Yeh R, Lee J, Nikhamin Y, Punzalan B, Punzalan B, La Perle K, Larson SM, Sadelain M, Brentjens RJ. Sensitive in vivo imaging of T cells using a membrane-bound *Gaussia princeps* luciferase. *Nat Med.* 2009 Mar;15(3):338-44. Epub 2009 Feb 15. PubMed PMID: 19219023; PubMed Central PMCID: PMC2837150.
- 30: Kim SB, Sato M, Tao H. Split *Gaussia* luciferase-based bioluminescence template for tracing protein dynamics in living cells. *Anal Chem.* 2009 Jan 1;81(1):67-74. PubMed PMID: 19061336.
- 31: Kim SB, Sato M, Tao H. Circularly permuted bioluminescent probes for illuminating ligand-activated protein dynamics. *Bioconjug Chem.* 2008 Dec;19(12):2480-6. PubMed PMID: 19049389.
- 32: Kim HJ, Chung JK, Hwang do W, Lee DS, Kim S. In vivo imaging of miR-221 biogenesis in papillary thyroid carcinoma. *Mol Imaging Biol.* 2009 Mar-Apr;11(2):71-8. Epub 2008 Nov 22. PubMed PMID: 19030936.
- 33: Kim HJ, Kim YH, Lee DS, Chung JK, Kim S. In vivo imaging of functional targeting of miR-221 in papillary thyroid carcinoma. *J Nucl Med.* 2008 Oct;49(10):1686-93. Epub 2008 Sep 15. PubMed PMID: 18794255.
- 34: Inouye S, Sahara Y. Soluble protein expression in *E. coli* cells using IgG-binding domain of protein A as a solubilizing partner in the cold induced system. *Biochem Biophys Res Commun.* 2008 Nov 21;376(3):448-53. Epub 2008 Sep 21. PubMed PMID: 18789309.
- 35: Takenaka Y, Masuda H, Yamaguchi A, Nishikawa S, Shigeri Y, Yoshida Y, Mizuno H. Two forms of secreted and thermostable luciferases from the marine copepod crustacean, *Metridia pacifica*. *Gene.* 2008 Dec 1;425(1-2):28-35. Epub 2008 Aug 7. PubMed PMID: 18725274.
- 36: Asokan A, Johnson JS, Li C, Samulski RJ. Bioluminescent virion shells: new tools for quantitation of AAV vector dynamics in cells and live animals. *Gene Ther.* 2008 Dec;15(24):1618-22. Epub 2008 Jul 31. PubMed PMID: 18668144; PubMed Central PMCID: PMC2684674.
- 37: Ketteler R, Seed B. Quantitation of autophagy by luciferase release assay. *Autophagy.* 2008 Aug 16;4(6):801-6. Epub 2008 Jun 3. PubMed PMID: 18641457; PubMed Central PMCID: PMC2910585.

- 38: Goerke AR, Loening AM, Gambhir SS, Swartz JR. Cell-free metabolic engineering promotes high-level production of bioactive *Gaussia princeps* luciferase. *Metab Eng*. 2008 May-Jul;10(3-4):187-200. Epub 2008 May 2. PubMed PMID: 18555198.
- 39: Ruecker O, Zillner K, Groebner-Ferreira R, Heitzer M. *Gaussia*-luciferase as a sensitive reporter gene for monitoring promoter activity in the nucleus of the green alga *Chlamydomonas reinhardtii*. *Mol Genet Genomics*. 2008 Aug;280(2):153-62. Epub 2008 May 31. PubMed PMID: 18516621.
- 40: Shao N, Bock R. A codon-optimized luciferase from *Gaussia princeps* facilitates the in vivo monitoring of gene expression in the model alga *Chlamydomonas reinhardtii*. *Curr Genet*. 2008 Jun;53(6):381-8. Epub 2008 Apr 12. PubMed PMID: 18408930; PubMed Central PMCID: PMC2413079.
- 41: Ketteler R, Sun Z, Kovacs KF, He WW, Seed B. A pathway sensor for genome-wide screens of intracellular proteolytic cleavage. *Genome Biol*. 2008 Apr 3;9(4):R64. PubMed PMID: 18387192; PubMed Central PMCID: PMC2643935.
- 42: Michelini E, Cevenini L, Mezzanotte L, Ablamsky D, Southworth T, Branchini BR, Roda A. Combining intracellular and secreted bioluminescent reporter proteins for multicolor cell-based assays. *Photochem Photobiol Sci*. 2008 Feb;7(2):212-7. Epub 2008 Jan 14. PubMed PMID: 18264589.
- 43: Hewett JW, Nery FC, Niland B, Ge P, Tan P, Hadwiger P, Tannous BA, Sah DW, Breakefield XO. siRNA knock-down of mutant torsinA restores processing through secretory pathway in DYT1 dystonia cells. *Hum Mol Genet*. 2008 May 15;17(10):1436-45. Epub 2008 Feb 7. PubMed PMID: 18258738; PubMed Central PMCID: PMC2861568.
- 44: Wurdinger T, Badr C, Pike L, de Kleine R, Weissleder R, Breakefield XO, Tannous BA. A secreted luciferase for ex vivo monitoring of in vivo processes. *Nat Methods*. 2008 Feb;5(2):171-3. Epub 2008 Jan 20. PubMed PMID: 18204457; PubMed Central PMCID: PMC2699561.
- 45: Lee JY, Kim S, Hwang do W, Jeong JM, Chung JK, Lee MC, Lee DS. Development of a dual-luciferase reporter system for in vivo visualization of MicroRNA biogenesis and posttranscriptional regulation. *J Nucl Med*. 2008 Feb;49(2):285-94. Epub 2008 Jan 16. PubMed PMID: 18199619.
- 46: Michelini E, Cevenini L, Mezzanotte L, Ablamsky D, Southworth T, Branchini B, Roda A. Spectral-resolved gene technology for multiplexed bioluminescence and high-content screening. *Anal Chem*. 2008 Jan 1;80(1):260-7. Epub 2007 Nov 22. PubMed PMID: 18031021.

- 47: Inouye S, Sahara Y. Identification of two catalytic domains in a luciferase secreted by the copepod *Gaussia princeps*. *Biochem Biophys Res Commun*. 2008 Jan 4;365(1):96-101. Epub 2007 Nov 5. PubMed PMID: 17981153.
- 48: Suzuki T, Usuda S, Ichinose H, Inouye S. Real-time bioluminescence imaging of a protein secretory pathway in living mammalian cells using *Gaussia luciferase*. *FEBS Lett*. 2007 Oct 2;581(24):4551-6. Epub 2007 Aug 27. PubMed PMID: 17825828.
- 49: Cheng G, Davis RE. An improved and secreted luciferase reporter for schistosomes. *Mol Biochem Parasitol*. 2007 Oct;155(2):167-71. Epub 2007 Jul 1. PubMed PMID: 17681388.
- 50: Badr CE, Hewett JW, Breakefield XO, Tannous BA. A highly sensitive assay for monitoring the secretory pathway and ER stress. *PLoS One*. 2007 Jun 27;2(6):e571. PubMed PMID: 17593970; PubMed Central PMCID: PMC1892804.
- 51: Venisnik KM, Olafsen T, Gambhir SS, Wu AM. Fusion of *Gaussia luciferase* to an engineered anti-carcinoembryonic antigen (CEA) antibody for in vivo optical imaging. *Mol Imaging Biol*. 2007 Sep-Oct;9(5):267-77. PubMed PMID: 17577599.
- 52: Morlighem JE, Petit C, Tzertzinis G. Determination of silencing potency of synthetic and RNase III-generated siRNA using a secreted luciferase assay. *Biotechniques*. 2007 May;42(5):599-600, 602, 604-6. PubMed PMID: 17515198.
- 53: Hewett JW, Tannous B, Niland BP, Nery FC, Zeng J, Li Y, Breakefield XO. Mutant torsinA interferes with protein processing through the secretory pathway in DYT1 dystonia cells. *Proc Natl Acad Sci U S A*. 2007 Apr 24;104(17):7271-6. Epub 2007 Apr 11. PubMed PMID: 17428918; PubMed Central PMCID: PMC1855419.
- 54: Knappskog S, Ravneberg H, Gjerdrum C, Trösse C, Stern B, Pryme IF. The level of synthesis and secretion of *Gaussia princeps luciferase* in transfected CHO cells is heavily dependent on the choice of signal peptide. *J Biotechnol*. 2007 Mar 10;128(4):705-15. Epub 2006 Dec 12. PubMed PMID: 17316861.
- 55: Trösse C, Ravneberg H, Stern B, Pryme IF. Vectors encoding seven oikosin signal peptides transfected into CHO cells differ greatly in mediating *Gaussia luciferase* and human endostatin production although mRNA levels are largely unaffected. *Gene Regul Syst Bio*. 2007 Dec 11;1:303-12. PubMed PMID: 19936096; PubMed Central PMCID: PMC2759145.
- 56: Remy I, Michnick SW. A highly sensitive protein-protein interaction assay based on *Gaussia luciferase*. *Nat Methods*. 2006 Dec;3(12):977-9. Epub 2006 Nov 12. PubMed PMID: 17099704.

57: Wiles S, Ferguson K, Stefanidou M, Young DB, Robertson BD. Alternative luciferase for monitoring bacterial cells under adverse conditions. *Appl Environ Microbiol.* 2005 Jul;71(7):3427-32. PubMed PMID: 16000745; PubMed Central PMCID: PMC1169068.

58: Tannous BA, Kim DE, Fernandez JL, Weissleder R, Breakefield XO. Codon-optimized *Gaussia* luciferase cDNA for mammalian gene expression in culture and in vivo. *Mol Ther.* 2005 Mar;11(3):435-43. PubMed PMID: 15727940.

59: Verhaegent M, Christopoulos TK. Recombinant *Gaussia* luciferase. Overexpression, purification, and analytical application of a bioluminescent reporter for DNA hybridization. *Anal Chem.* 2002 Sep 1;74(17):4378-85. PubMed PMID: 12236345.